

Polyclonal Antibody against Human CRP

Catalog Number: 11120 Size: 100 µg Host: Rabbit

Introduction to the Molecule

C-reactive protein (CRP) is a circulating protein mainly secreted from the liver. This acute phase protein consists of five identical non-glycosylated subunits of 23 kDa, that give rise to a symmetrically arranged globular protein with molecular weight of approximately 120 kDa¹. It has long been recognized that CRP is closely related to immunology, inflammation and host defence, as a result it has been used as an inflammatory marker. There is accumulating evidence suggesting the important role that CRP plays in mediating cardiovascular diseases (CVD) and type 2 diabetes²⁻⁴. Normally, CRP is presenting only in a trace amount in circulation (<1 μ g/mL)^{5,6} but can increase over 1,000-fold under acute inflammatory state. Individual with blood CRP levels <1 μ g/mL, 1-3 μ g/mL and >3 μ g/mL is considered to have low, moderate and high risk, respectively, of CVD and myocardial infraction⁷. Therefore, blood CRP level has become a promising measure of CVD risk^{8,9}.

Purification

Rabbit crude IgG was purified by protein-G column.

Immunogen

Recombinant Human C-reactive Protein (CRP) expressed in E.coli.

Specificity

The antibody detects human CRP.

Formulation & Storage

Liquid in phosphate-buffered saline (PBS). Store at -20°C for less than one week. For long-term storage, aliquot and freeze at -70°C. Avoid repeated freeze/defrost cycles.

Application/Usage

ELISA - This antibody can be used at 0.1-0.2 μ g/mL with the appropriate secondary reagents to detect human CRP.

Quality Control Test

BCA to determine quantity of the antibody.

References

[1] Thompson D., Pepys M.B. and Wood S.P. (1999) Structure, 7, 169-177.

[2] Festa A, D'Agostino R. Jr., Tracy R.P. and Haffner S.M. (2002) Diabetes, 51, 1131-1137.

- [3] Verma S. and Yeh E.T. (2003) Am J Physiol, 285, R1253-R1258.
- [4] Jialal I., Devaraj S. and Venugopal S.K. (2004) Hypertension, 44, 6-11.
- [5]Kindmark C.O. (1972) Scand J Clin Lab Invest, 29, 407-411.
- [6] Macy E.M., Hayes T.E. and Tracy R.P. (1997) Clin Chem, 43, 52-58.
- [7]Ridker P.M. (2004) Am Heart Hosp J, 2 (4 Suppl 1), 4-9.
- [8]Benzaquen I.R., Yu H. and Rifai N. (2002) Crit Rev Clin Lab Sci, 39, 459-497.
- [9]Pearson T.A. et al., (2003) Circulation, 107, 499-511.

Website: www.immunodiagnostics.com.hkE-mail: info@immunodiagnostics.com.hk (HK) / info@immunodiagnostics.ca (Canada)Tel: +852 3502 2780 (HK) / +1-437-886-5136 (Canada)